Memorandum

Date: October 16, 2017

To: Kris Cano, Forensic Laboratory Manger *WCW 1/5/18*
Mindy Raines, Forensic Scientist Supervisor *MR*

From: Allan Kosecki, Blood Alcohol Technical Leader

Re: Incorrect Target Value for Control ACQ Science 407041529/3

The target value supplied on the certificate of analysis (COA) for the ACQ Science blood control lot 407041529/3 is 0.198 g/dL. The target value for the previous lot, 402101060/6, was 0.200 g/dL. When the new lot was first used in case work, run 10Feb17, the analyst did not update the value to 0.198 g/dL. This error was carried forward in all runs from February 10, 2017 through October 3, 2017. As a result, all data generated for the blood alcohol runs completed during this period have the incorrect target value for control ACQ Science 407041529/3.

With a target value of 0.198 g/dl, the acceptable range for agreement within five percent is less than 0.208 g/dL and greater than 0.188 g/dL. All analyses of the whole blood control lot 407041529/3 in the runs conducted between February 10, 2017 and October 3, 2017 were within this range.

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EtOH WH 2.0 g/l – For in vitro diagnostic use

Ethanol control in whole blood

Application

This material should be used in accordance with the laboratory's operating procedures for instrument calibration or as a control material.

User guide

This ACQ Science EtOH WH requires no additional preparation and is ready for use.

Assigned value

The assigned ethanol concentration was determined by 3 independent laboratories, each accredited to DIN EN 17025. Repeat determinations were carried out daily on 5 days using Gas Chromatography.

Storage and stability

Storage: 2° to 8° C

Stability:

- Sealed container, stored in the dark: see expiration date on the package.
- Stored in the dark tightly capped: see expiration date on package

Precautions

All materials of human origin should be considered as potentially infectious and treated with the same care as patient specimens. Each individual original blood unit used for the production of the control was tested for the following antigens and antibodies: HBsAG, anti-HIV-1, anti-HIV-2, anti-HBc and anti-HCV and found to be negative.

Lot: 402101060
 Order no.: WH20-015 (10 x 1,5 ml)
 WH20-030 (10 x 3,0 ml)
 Version: 3 – 201303

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GI_EtOHWH_20_402101060_De_En_20130325.doc

¹ Confidence ranges – measured values

The confidence interval indicates the range in which the target value is located with a significance level of 95%.

² Confidence ranges – German forensic directives

[EtOH] ≥ 1,0 g/l → ± 5% from the target value
 [EtOH] < 1,0 g/l → ± 0.05 g/l from the target value

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³ Confidence ranges – Directive of the German Medical Association

0,2 < [EtOH] ≤ 0,6 g/l → ± 15 % from the target value
 0,6 < [EtOH] ≤ 5,0 g/l → ± 9 % from the target value

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EtOH WH 2,0 g/l – For in vitro diagnostic use

Ethanol control in whole blood

Application

This material should be used in accordance with the laboratory's operating procedures for instrument calibration or as a control material.

User guide

This ACQ Science EtOH WH requires no additional preparation and is ready for use.

Assigned value

The assigned ethanol concentration was determined by 3 independent laboratories each accredited to DIN EN 17025. Repeat determinations were carried out daily on 5 days using Gas Chromatography.

Storage and stability

Storage: 2° to 8° C

Stability:

- Sealed container, stored in the dark: see expiration date on the package.
- Stored in the dark tightly capped: see expiration date on package

Precautions

All materials of human origin should be considered as potentially infectious and treated with the same care as patient specimens. Each individual original blood unit used for the production of the control was tested for the following antigens and antibodies: HBsAG, anti-HIV-1, anti-HIV-2, anti-HBc and anti-HCV and found to be negative.

Lot: 407041529
 Order no.: WH20-015 (10 x 1,5 ml)
 WH20-115 (100 x 1,5 ml)
 WH20-030 (10 x 3,0 ml)
 Version: 2 - 201510

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GI_EtOHWH_20_407041529_20151029

¹ Confidence ranges – measured values

The confidence interval indicates the range in which the target value is located with a significance level of 95%.

² Confidence ranges – German forensic directives

[EtOH] ≥ 1,0 g/l → ± 5% from the target value
 [EtOH] < 1,0 g/l → ± 0,05 g/l from the target value

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ACQ Science GmbH